





TRAIL, IP-10, CRP Host-protein Signature score distinguishes between Viral and Bacterial Infection in Sepsis Patients.

BSAC Spring Meeting, May 11th

Jeroen Stas, Medical Affairs Manager, MeMed Diagnostics.




Investigators of Apollo (sub-)study

- Sergey M. Motov, M.D. Maimonides Medical Center, Emergency Medicine. Brooklyn, US.
- Cesar A. Arias, M.D., Ph.D. University of Texas Health Science Center (UTHealth). Houston, US.
- Richard G. Bachur, M.D. Boston Children's Hospital. Boston, US.
- Susanna Esposito, M.D. Pediatric Clinic, Pietro Barilla Children's Hospital. Parma, Italy.
- Salim Halabi, M.D. Carmel Medical Center. Haifa, Israel.
- Sheldon L. Kaplan, M.D. Texas Children's Hospital. Houston, US.
- Adi Klein, M.D. Hillel Yaffe Medical Center. Hadera, Israel.
- Richard Rothman, M.D., Ph.D. Johns Hopkins University. Baltimore, US.
- Leticia M. Ryan, M.D., M.P.H. Johns Hopkins University. Baltimore, US.
- Shachaf Shiber, M.D. Rabin Medical Center. Petah Tikva, Israel.
- Tobias Tenenbaum, M.D. Sana Klinikum Lichtenberg. Berlin, Germany.
- Alexandra Weissman, M.D. University of Pittsburgh Medical Center. Pittsburgh, US.

DISCLAIMER: MeMed funded the Apollo study with payments to participating institutions. I am a MeMed employee.

Changing Clinical Diagnostics by measuring the host-response




What Bug?




Inaccessible infection sites



False Alarms due to colonization



Often, no pathogens are detected



Poor performance in emerging pathogens



How is the patient responding?



Measure & quantify relevant immune response biomarkers.




Integration of biomarkers in a clinical score instead of using cutoffs



Small blood volume (100µl serum) in 15 minutes

Background: Host-protein Signature score.



TRAIL




CRP

IP-10


TNF-Related Apoptosis-Inducing Ligand

C-Reactive Protein

Interferon Gamma-Induced Protein 10



Background: Host-protein Signature score.



Validation: Host-protein Signature score

Study Name	Patients #	Clinical Setting	Age	Clinical Syndrome	BV Performance
Curiosity / ISR [PLOS One, 2015]	1,002	ED, wards	Pediatrics and Adults, ≥3mo	Suspicion of acute infection	AUC 93%
Pathfinder* / ISR, SWISS [Pediatrics, 2017]	597	ED, wards	Pediatrics, ≥3mo, ≤18yr	Pneumonia, FWS, blood drawn for serology	AUC 92%
Opportunity / ISR,NL [The Lancet, 2016]	777	ED, wards	Pediatrics, 2-60 months	Symptoms of LRTI or FWS	AUC 90%
Rosetta / ISR [Data on File]	550	ED	Pediatrics, ≥3mo, ≤18yr	Symptoms of RTI or FWS	AUC 91%
Autopilot / GER, ITA [CMI, 2021]	1,140	ED	Pediatrics, ≥3mo, ≤18yr	Suspicion of RTI or FWS	AUC 94%
Observer / ISR [Data on File]	583	ED	Adults, >18 years	Suspicion of LRTI or pneumonia	AUC 98%
Apollo / US, ISR, GER*, ITA* [Data on File]	1,016	ED, UCC	Pediatrics and Adults, > 3mo	Suspicion of acute bacterial or viral infection	AUC 97%

*Retrospective

The Apollo study



Regulatory

Work in US patient population and design new FDA regulatory pathway for host-response testing.



Validate

MeMed BV test together with the MeMed Key platform.



Compare

Diagnostic performance to other inflammatory biomarkers and their combinations.


Apollo study

Inclusion:

- Older than 90 days of age
- Suspected viral or bacterial infection
- Febrile within last 7 days
- Duration of illness less than 7 days

Exclusion:

- Unrelated febrile episode within 2 weeks
- Immunocompromised



Clinical setting

Emergency department and / or urgent care centers (no hospital-admitted inpatients).



Expert panel adjudication

Using the full patient case file + 28-day follow-up data.



The Apollo study: sub-analysis on Sepsis patients

- Sub-study objective:
 - “Evaluate the MeMed BV’s ability to differentiate viral from bacterial infection in Sepsis patients”
- Sepsis definition:
 - “Two or more SIRS criteria”.
 - Still commonly used in the US, despite Sepsis-3 definition.
 - Expert adjudication: Viral or bacterial classification by at least 2/3 experts to assign same etiology label with confidence $\geq 90\%$ or all 3 assign with confidence $\geq 70\%$.

SIRS criteria for children:

Age Group	Pediatric SIRS Criteria (≥ 1 of the criteria from Column 1 AND Column 2)							Systolic Blood Pressure (mmHg)	
	Column 1 (≥ 1 of the below criteria)				Column 2 (≥ 1 of the below criteria)				Cardiovascular Dysfunction
	Core Temperature ($^{\circ}\text{C}$)		Leukocyte Count (Leukocytes $\times 10^3/\text{mm}^3$) ³		Heart Rate (Beats/Min) ¹				
	Hypothermia	Hyperthermia	Leukopenia	Leukocytosis	Bradycardia	Tachycardia	Respiratory Rate ² (Breaths/Min)		
0 days to 1 wk	<36	>38.5	NA	>34	<100	>180	>50	<65	
1 wk to 1 mo	<36	>38.5	<6	>19.5	<100	>180	>40	<75	
1 mo to 1 yr	<36	>38.5	<6	>17.5	<90	>180	>34	<100	
2-5 yrs	<36	>38.5	<6	>15.5	NA	>140	>22	<94	
6-12 yrs	<36	>38.5	<4.5	>13.5	NA	>130	>18	<104	
13 to <18 yrs	<36	>38.5	<4.5	>11	NA	>110	>14	<117	

SIRS criteria for adults:

Box 1. SIRS (Systemic Inflammatory Response Syndrome)

Two or more of:

Temperature $>38^{\circ}\text{C}$ or $<36^{\circ}\text{C}$


Heart rate $>90/\text{min}$

Respiratory rate $>20/\text{min}$ or $\text{PaCO}_2 <32 \text{ mm Hg (4.3 kPa)}$

White blood cell count $>12\,000/\text{mm}^3$ or $<4000/\text{mm}^3$ or $>10\%$ immature bands

From Bone et al.⁹


Sub-analysis: Results – Patient Flow




**This was a 9-year-old child presenting within 2 days of fever with a sore throat and pharyngeal erythema; rapid strep A test was negative, antibiotic was prescribed, and the child was discharged.

HR = heart rate; RR = respiratory rate

Sub-analysis: Results – Patient Characteristics



Median age of children was 2.4 years (IQR: 1.4-5.4); adults was 41.8 years (IQR: 29.2-61)



Median duration of hospitalization was 4 days (IQR: 3-6)

Sub-analysis: Results – Clinical Syndromes

URTI discharge diagnoses included:

Acute Otitis Media; Tonsillitis; Upper Respiratory Infections; Peritonsillar Abscess; Pharyngitis; Sinusitis; Aphthous Stomatitis; Coronavirus (Covid19) Infection; Scarlet Fever; Flu/Influenza; Flu Like Illness/Symptoms; Herpangina; Laryngitis; Respiratory Viral Syndrome; Tracheitis

LRTI discharge diagnoses included:

Bronchiolitis; Acute Bronchitis; Acute Respiratory Failure; Bronchopneumonia; COPD Exacerbation; Community Acquired Pneumonia; Lobar Pneumonia; Pneumonia

UTI discharge diagnoses included:


Acute Cystitis; Pyelonephritis; Urinary Tract Infection

CNS discharge diagnoses included:

Meningitis


The 'Other' category included:

Abdominal Pain; Abscess; Appendicitis; Asthma; Cellulitis; Febrile Convulsions; Fever; Gastroenteritis; Headache; Unspecified viral infection



Patients could be included in more than one clinical syndrome except for Other

Sub-analysis: MeMed BV performance



MeMed BV Score Cut-offs:

- < 35 - viral infection
- $35 \leq \text{score} \leq 65$ – equiv. (7.9%)
- > 65 - bacterial infection

Procalcitonin Cut-offs

Sepsis algorithm:

- $\geq 0.5\text{ng/mL}$ - bacterial infection
- $< 0.5\text{ng/mL}$ - viral infection

LRTI algorithm:

- $\geq 0.25\text{ng/mL}$ - bacterial infection
- $< 0.25\text{ng/mL}$ - viral infection

(nm)

-70.5)

-77.2)

-87.8)

Apollo Sepsis sub-analysis - Conclusions

- MeMed BV accurately distinguished viral from bacterial infection in sepsis patients.
- This new triage tool has potential to help with timely identification of bacterial infection, enabling prompt treatment.
- MeMed BV accurately rules out bacterial infection, allowing antibiotic overuse to be minimized.



Thank you

